

KAPI-GEL®

tracheostoma spacer

- ✔ Ideal for difficult stoma anatomy. Helps to 'fill in' the stoma area, allowing the comfortable use of Laryngectomy tubes, stoma studs, Heat Moisture exchange/filter holders.
- ✔ Acts as a cushion behind the above style of tubes, particularly useful if one is digitally occluding the stoma.
- ✔ Due to its transparency, Kapi-Gel is not noticeable around the stoma
- ✔ Will help with the location of stoma studs and laryngectomy tubes
- ✔ Promotes skin care around the stoma and helps reduce soreness around the stoma and tube site.
- ✔ May help reduce air leakage around the tube.
- ✔ Kapi-Gel is a tacky laminated silicone gel available in two thicknesses 3mm and 5mm and in two hole dimensions 8mm and 12mm. Spacer diameter 43mm.
- ✔ Available in packs of 10.



NEW!



Easy to remove release paper



Position over stoma*



Kapi-Gel spacer sitting behind Trachi-Naze Plus stud



Kapi-Gel spacer on the shaft of a Trachi-Naze Plus laryngectomy tube



Kapi-Gel spacer on the shaft of a Trachi-Naze Plus stoma stud

ORDER REFERENCES

LA TNG 2008	Kapi-Gel THIN 3mm	8mm pack of 10
LA TNG 2012	Kapi-Gel THIN 3mm	12mm pack of 10
LA TNG 4008	Kapi-Gel THICK 5mm	8mm pack of 10
LA TNG 4012	Kapi-Gel THICK 5mm	12mm pack of 10



ISO 9000 / 2000
ISO 13485
Cert. 257



Kapitex Healthcare Ltd

Kapitex House, 1 Sandbeck Way, Wetherby, West Yorkshire LS22 7GH England

Tel: +44 (0)1937 580211 Fax: +44 (0)1937 580796 Email: sales@kapitex.com

www.kapitex.com

*It may be easier to place Kapi-Gel on the tube/stud shaft before placing around stoma



KAPITEX®

The Specialists in Laryngectomy,
Tracheostomy, Dysphagia and
Airway Management